



Perfect **Skn** Lite



This is a fantastic beauty system that treats ageing, vascular and acne skins all rolled up into one.

An advanced beauty system that uses LED photobiology that emits light which acts on skin to raise cellular activity and improve skin metabolism.



Why Perfect Skn Lite For You?

- Easy to use system
- Promises fantastic results ... and gives them
- Versatile system to treat ageing, vascular and acne skin
- Easy to sell facial with high profitability
- All light waves incorporated into one single machine

Perfect Skn Lite

Red Light - The red light is absorbed into skin cells which increases the energy of the cell, accelerating normal cellular process and making the skin cells work more efficiently. This has the effect of stimulating the collagen thus causing skin rejuvenation. This treatment also effects whitening, speckle and vascular removing, reducing acne and recovering injured skin.

Blue Light - The blue light targets and neutralises acne bacteria thus eliminating acne. Acne treatment would be carried out by alternating Red Light and Blue Light twice a week for four weeks. The Red Light stimulates anti-inflammatory properties to reduce redness, promote skin healing whilst the Blue Light eliminates acne bacteria.

Mixed Light - The mix of blue and red light together stimulates collagen and elastin production, eliminates fine wrinkles, enables skin to look tighter, more even toned and clearer.

Perfect Skn Lite



specifications		
Light Source	Output Wavelength	Spot Size
LED photobiology light	Red: 640nm±5nm Blue: 415nm±5nm	350 x 240mm
Surface Temperature of Treatment Head	Output Volt	Input Power
≤40°C	80mw/cm	220V 50Hz 110V 60Hz
Size	Weight	
630x245x560mm	12.5kg	